



Sample Description

Analyzed date: 2016/5/17

Analyst: Shen

Sample: P-16

M.W.: 1658.78

Lot. No.: P160509-LL654371

Instrument

SHIMADZU LCMS-2020

Probe:

ESI

Probe Bias:

+ 4.5kv

Nebulizer Gas Flow:

1.5L/min

Detector:

1.2kv

CDL:

-20.0v

T. Flow:

0.2ml/min

CDL Temp:

250 °C

B. Conc.:

50% H2O/50% ACN

Block Temp:

400 °C