


HPLC REPORT

Sample: G-16
 Analyst: YW
 Lot. No.: P150517-SY654820
 Column: 4.6×250mm, Kromasil 100-5-C18
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	15%	85%
25.0min	40%	60%
25.1min	100%	0%
30.0min	Stop	

Volume: 5µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min

Analyzed date: 2015-05-21



Peak	Time	Height	Area	Conc.
1	7.295	763.581	7091.609	0.4221
2	13.787	3164.168	33810.320	2.0127
3	14.358	9009.279	76477.328	4.5525
4	14.608	106556.813	1460221.625	86.9236
5	15.437	3056.385	53483.613	3.1838
6	15.725	2129.402	23174.820	1.3795
7	15.977	1894.928	22459.006	1.3369
8	16.485	400.206	3171.903	0.1888
Total				100.000