


HPLC REPORT

Sample: L-14 Analyzed date: 2015-8-23
Analyst: LY
Lot. No.: P150815-YS121564
Column: Kromasil-C18, 4.6*250mm, 5µm
Solvent A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

	A	B
0.0min	16%	84%
25.0min	41%	59%
25.1min	100%	0%
30.0min	Stop	

Volume: 10µl
Wavelength: 214nm
Flow rate: 1.0ml/min



Peak No	Ret Time	Height	Area	Conc.
1	10.548	1451.985	16834.420	0.5463
2	10.740	3828.115	13283.243	0.4311
3	10.905	304561.750	2936344.000	95.2878
4	10.905	10043.666	115091.484	3.7349
Total				100.00