


HPLC REPORT

Sample: S-16
 Analyst: LLR
 Lot. No.: P140725-HX599248
 Column: Kromasil C18, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	10%	90%
25.0min	35%	65%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min

Analyzed date: 2014-08-02



Rank	Time	Area	Conc.
1	10.812	109293.344	9.0953
2	10.997	1045441.875	87.0004
3	10.997	46917.004	3.9044

Total 1201652.223 100