


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2015-8-31
 Analyst: SHAO
 Sample: T-15
 M.W.: 1836.23
 Lot. No.: P150824-SY673157

Instrument

SHIMADZU LCMS-2020
 Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp.: 250°C
 Block Temp.: 200°C
 Probe Bias: + 4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. Conc.: 50%H2O/50%ACN