

HPLC REPORT

Sample: C-16

Analyzed date: 2015-08-24

Analyst: ZK

Lot. No.: P150821-SY673183

Column: Kromasil 100-5C18, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water


Gradient :

	A	B
0.0min	25%	75%
25.0min	50%	50%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Quantity	Area	Height
1	12.350	8.1085	407522	7218
2	12.818	89.1524	4480705	421173
3	13.238	2.7391	137663	5144
Total		100	5025890	433535