


HPLC REPORT

Sample: L-14
 Analyst: YW
 Lot. No.: P150515-SY654967
 Column: 4.6×250mm, Kromasil 100-5-C18
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	22%	78%
25.0min	47%	53%
25.1min	100%	0%
30.0min	Stop	

Volume: 5µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min

Analyzed date: 2015-05-21



Peak	Time	Height	Area	Conc.
1	9.013	5720.375	100968.250	0.6047
2	10.003	24146.123	705601.688	4.2259
3	10.523	38951.699	1195047.625	7.1573
4	11.317	742813.625	14648583.000	87.7321
5	11.317	4572.983	46750.688	0.2800

Total				100.000
-------	--	--	--	---------