


# HPLC REPORT

Sample: C-11 Analyzed date: 2016-05-22  
Analyst: ZK  
Lot. No.: P160514-SY655831  
Column: Kromasil 100-5C18, 4.6\*250mm, 5µm  
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
Gradient :

|         | A    | B    |
|---------|------|------|
| 0.0min  | 15%  | 85%  |
| 25.0min | 40%  | 60%  |
| 25.1min | 100% | 0%   |
| 30.0min |      | Stop |

Volume: 20µl  
Wavelength: 220nm  
Flow rate: 1.0ml/min



| Rank  | Time   | Quantity | Area    | Height |
|-------|--------|----------|---------|--------|
| 1     | 12.636 | 14.6953  | 929552  | 122550 |
| 2     | 12.922 | 85.3047  | 5395956 | 678791 |
| Total |        | 100      | 6325508 | 801341 |