


MASS SPECTROMETRY REPORT

19-Apr-2018
10:36:07

Probe: ESI
Cone: 50v
Desovation Temp :350

Capillary:3.00KV
Extractor: 5v
Gas Flow : 350
Scan ES+
6.34e7

180419-9 1105 (20.424)



Sample Description

Analyzed date: 2016-04-19
Analyst: YU
Sample: V-12
M.W.: 1027.21
Lot. No.: P160413-CJ649278

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50% H2O/50% ACN