

CAPP

Accessories

MECHANICAL PIPETTES

CAPP offers a comprehensive selection of practical accessories to complement its full line of mechanical pipettes. Uniquely designed and quality made they combine design, style and simplicity to deliver an optimal storage and safety solutions for modern laboratories.

Precalibrated interchangeable adjustable and fixed volume knobs set **CAPP** pipettes aside from any other liquid handling pipette on a market by allowing users to use one pipette with many fixed volume knobs to substantially reduce cost.

Now, instead of four fixed volume pipettes, for example, you can buy one **CAPP** pipette and four inexpensive knobs to achieve the same precise results.

- **Stackable CAPP Combi Stand** can be connected to one another to create infinite support from one pipette to as many as user might need. C-01 comes with 2 modules for support of up to 3 mechanical pipettes
- **CAPP Carousel S-06** provides ultimate space saving solution for up to 6 mechanical pipettes
- **T-01 CAPP Timer** is a simple, hence reliable tool with audible alarm feature
- **Shaft filters** provide an economical solution to eliminate “splash effect” and will substantially reduce risk of sample and shaft contamination



Ordering Information

Pipette Filters

Cat. No.	Description
C1000-1-FT	Disposable filter for use with CAPP ecopipette 100-1000µL, bag w/ 50 pcs.
C5000-1-FT	Disposable filter for use with CAPP ecopipette 1-5mL, bag w/ 25 pcs.
B510ML-1-FT	Disposable filter for use with CAPPBravo pipettes, 1-5mL & 1-10mL, bag w/ 25 pcs.