


## Ordering information:

| Code    | Description                                     | Unit      | Packaging   |
|---------|---|-----------|-------------|
| CA40505 | 30 ml reservoirs for 8- and 16-channel pipettes | Pack/50   | Non-Sterile |
| CA40506 | 30 ml reservoirs for 8- and 16-channel pipettes | Pack/5x10 | Sterile     |
| CA40510 | 40 ml reservoirs for 12-channel pipettes        | Pack/50   | Non-Sterile |
| CA40511 | 40 ml reservoirs for 12-channel pipettes        | Pack/5x10 | Sterile     |



## CappAero™ Multi

| Standard | Type       | Volume    | Colour   |
|----------|------------|-----------|----------|
| C10-8    | 8-channel  | 0.5-10 µl | ○ White  |
| C20-8    | 8-channel  | 2-20 µl   | ● Grey   |
| C50-8    | 8-channel  | 5-50 µl   | ● Red    |
| C100-8   | 8-channel  | 10-100 µl | ● Orange |
| C200-8   | 8-channel  | 20-200 µl | ● Yellow |
| C300-8   | 8-channel  | 30-300 µl | ● Black  |
| C10-12   | 12-channel | 0.5-10 µl | ○ White  |
| C20-12   | 12-channel | 2-20 µl   | ● Grey   |
| C50-12   | 12-channel | 5-50 µl   | ● Red    |
| C100-12  | 12-channel | 10-100 µl | ● Orange |
| C200-12  | 12-channel | 20-200 µl | ● Yellow |
| C300-12  | 12-channel | 30-300 µl | ● Black  |
| C02-16   | 16-channel | 0.2-2 µl  | ● Violet |
| C10-16   | 16-channel | 0.5-10 µl | ○ White  |
| C50-16   | 16-channel | 5-50 µl   | ● Red    |



# CAPP®

## Capp A/S

Landbrugsvej 10 • DK 5260 Odense S  
Denmark  
tel. +45 6613 6140 • fax. +45 6613 2770  
[www.capp.dk](http://www.capp.dk) • [info@capp.dk](mailto:info@capp.dk)

Distributor