


CAPP *Origami*

REAGENT RESERVOIRS



Use your reagent most efficiently with **CAPP *Origami*** Reagent Reservoirs.

The **CAPP *Origami*** represents a new way of thinking when it comes to reagent reservoirs.

The stiff architecture procure a consistently sharp V-shape reservoir. This makes pipetting easier and assures an efficient use of the reagent.

They are packed, shipped and stored as flat sheets.

To use the reservoir, simply flick it up by pushing the sides together. This environmentally and storage friendly product will save bench space and money.

Ordering Information

Reagent Reservoirs: **CAPP *Origami***™ 30mL and 40mL reservoirs for multichannel pipetting

Cat. No.	Description
CA40505	CAPP Origami 30mL (8- and 16-channel pipettes), bag w/ 50 pcs.
CA40506	CAPP Origami 30mL (8- and 16-channel pipettes), Pre-sterile, 10 bags w/ 5 pcs. each
CA40510	CAPP Origami 40mL (12-channel pipettes), bag w/ 50 pcs.
CA40511	CAPP Origami 40mL (12-channel pipettes), Pre-sterile, 10 bags w/ 5 pcs. each

CAPP *Aero 96* Multichannel pipettes

Cat. No.	Description	Inaccuracy %	Imprecision %	Colour
C10-8	CAPP multi pipettes, 8-channel, 0.5-10µL	5.00 /2.00	3.60 /1.00	○
C20-8	CAPP multi pipettes, 8-channel, 2-20µL	5.00 /2.00	3.60 /1.00	●
C50-8	CAPP multi pipettes, 8-channel, 5-50µL	3.20 /1.40	1.60 /0.60	●
C100-8	CAPP multi pipettes, 8-channel, 10-100µL	2.50 /1.30	1.20 /0.50	●
C200-8	CAPP multi pipettes, 8-channel, 20-200µL	1.80 /1.20	0.80 /0.40	●
C300-8	CAPP multi pipettes, 8-channel, 30-300µL	1.66 /1.40	0.60 /0.30	●
C10-12	CAPP multi pipettes, 12-channel, 0.5-10µL	5.00 /2.00	3.60 /1.00	○
C20-12	CAPP multi pipettes, 12-channel, 2-20µL	5.00 /2.00	3.60 /1.00	●
C50-12	CAPP multi pipettes, 12-channel, 5-50µL	3.20 /1.40	1.60 /0.60	●
C100-12	CAPP multi pipettes, 12-channel, 10-100µL	2.50 /1.30	1.20 /0.50	●
C200-12	CAPP multi pipettes, 12-channel, 20-200µL	1.80 /1.20	0.80 /0.40	●
C300-12	CAPP multi pipettes, 12-channel, 30-300µL	1.66 /1.40	0.60 /0.30	●

