

Ordering Information



CAPP *Bravo* Starter Packs

Cat. No.	Description
SK-03-BB	Three pipettes free of choice, C-05 stand, 3 racks of tips and 1 timer

Pipette Stands

Cat. No.	Description
C-05	Pipette stand for up to 5 pipettes

CAPP *Bravo* Variable Volume

Cat. No.	Description	Inaccuracy %	Imprecision %	Colour
B02-1	CAPP pipette, variable vol. 0.1-2 μ L, Bravo	12.00 / 2.00	6.00 / 1.00	●
B10-1	CAPP pipette, variable vol. 0.5-10 μ L, Bravo	2.50 / 1.00	1.80 / 0.50	○
B20-1	CAPP pipette, variable vol. 2-20 μ L, Bravo	5.00 / 1.00	1.50 / 0.30	●
B50-1	CAPP pipette, variable vol. 5-50 μ L, Bravo	2.00 / 1.00	1.20 / 0.30	●
B100-1	CAPP pipette, variable vol. 10-100 μ L, Bravo	2.50 / 0.70	1.00 / 0.20	●
B200-1	CAPP pipette, variable vol. 20-200 μ L, Bravo	2.50 / 0.60	0.70 / 0.20	●
B1000-1	CAPP pipette, variable vol. 100-1000 μ L, Bravo	2.00 / 0.80	0.60 / 0.20	●
B5000-1	CAPP pipette, variable vol. 1-5mL, Bravo	1.20 / 0.60	0.30 / 0.15	●
B10000-1	CAPP pipette, variable vol. 1-10mL, Bravo	3.00 / 0.60	0.60 / 0.15	●

CAPP *Bravo* Fixed Volume (customized volumes available on request)

Cat. No.	Description	Inaccuracy %	Imprecision %
B01-1F	CAPP pipette, fixed vol. 1 μ L, Bravo	2.50	1.80
B02-1F	CAPP pipette, fixed vol. 2 μ L, Bravo	2.00	1.00
B05-1F	CAPP pipette, fixed vol. 5 μ L, Bravo	1.20	0.80
B10-1F	CAPP pipette, fixed vol. 10 μ L, Bravo	1.00	0.50
B20-1F	CAPP pipette, fixed vol. 20 μ L, Bravo	1.00	0.30
B50-1F	CAPP pipette, fixed vol. 50 μ L, Bravo	1.00	0.30
B100-1F	CAPP pipette, fixed vol. 100 μ L, Bravo	0.70	0.20
B200-1F	CAPP pipette, fixed vol. 200 μ L, Bravo	0.60	0.20
B500-1F	CAPP pipette, fixed vol. 500 μ L, Bravo	0.80	0.15
B1000-1F	CAPP pipette, fixed vol. 1000 μ L, Bravo	0.80	0.20
B5000-1F	CAPP pipette, fixed vol. 5mL, Bravo	0.60	0.15
B10000-1F	CAPP pipette, fixed vol. 10mL, Bravo	0.60	0.15