



EXPell pipette tips reload system is available for 10µL standard and extended length as well as for the 200µL and 1000-1250µL tips. Depending on the tip volume, it consists of towers with 5 or 10 insert plates containing 96 tips each.

The system is extremely easy to use: take off the tower cover and pick up the insert plate with tips. Place it into the empty rack. The rack has been refilled with 96 tips.

EXPell reload system enables multiple savings. It reduces the shipping cost significantly by cutting off up to 75% of the volume being shipped. The tips in an **EXPell** reload system are packed more compact compared to standard racked or bulk tips.

In the **EXPell** reload system the tips from the upper insert plate are partially placed into the tips of the lower plate, resulting in a slim tower, taking much less space than a single box with 10 racked tips.

CAPP EXPell reload system reduces the plastic waste. Each rack can be reused multiple times and apart from the insert plates with tips, there is no further plastic component in the **EXPell** reload system – only the recyclable carton packaging.

Further saving with **EXPell** reload system concerns the space. Thanks to the minimal footprint, the towers of the **EXPell** reload system take as less space on the benchtop and in the storage room as possible.

The system is compatible with all **CAPP** racks of corresponding tip volumes.

Spare, empty racks are available as accessories.

Ordering Information

Cat. No.	Description
5030035C	Expell 10µL, reload, 10x10x96 pcs.
5030034C	Empty rack for 10µL reload tips, 50 pcs.
5030055C	Expell 10µL XL, reload, 10x10x96 pcs
5030054C	Empty rack for 10µL XL tips, 50 pcs.
5030095C	Expell 200µL, reload, 10x10x96 pcs.
5030094C	Empty rack for 200µL reload tips, 50 pcs.
5130155C	Expell 1000µL (1250µL), reload, 10x5x96 pcs.
5130125C	Empty rack for 1000µL (1250µL) reload tips, 32 pcs.