

# EX Pell and EX PellPLUS

## COMPATIBILITY CHART

For tips' compatibility with other pipette brands, please contact [info@capp.dk](mailto:info@capp.dk)



5130070C 5030070C 5030072C 4130075C 5030075C 5030080C 5030095C	5030089C 5030090C 5130090C 5131090C	5130100C 5130105C 5130110C	5030114C 5030115C	5130130C 5130132C 4130135C 5130135C 5130140C 5130155C	5130149C 5130150C 5131150C	5130120C 5130121C 5130122C	5130119C 5130123C	5130160C 5130165C	5130190C 5130195C
2-20µL 5-50µL 10-100µL 20-200µL	2-20µL 5-50µL 10-100µL 20-200µL	2-20µL 5-50µL 10-100µL 20-200µL	2-20µL 5-50µL 10-100µL 20-200µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	1-5mL	1-10mL
5-50µL 10-100µL 20-200µL	5-50µL 10-100µL 20-200µL	5-50µL 10-100µL 20-200µL 30-300µL	5-50µL 10-100µL 20-200µL 30-300µL	100-1200µL	100-1200µL	100-1200µL	100-1200µL	500-5000µL	x
2-20µL 10-100µL 20-200µL	10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	x	1-10mL
2-20µL 10-100µL 20-200µL	20-200µL	20-200µL 15-300µL	20-200µL 15-300µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	x	x
2-20µL 10-100µL 20-200µL	20-200µL	10-100µL 20-200µL 30-300µL	10-100µL 20-200µL 30-300µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	x	1-10mL
2-20µL 5-50µL 20-200µL	5-50µL 20-200µL	5-50µL 20-200µL 30-300µL	5-50µL 20-200µL 30-300µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	x	x
2-20µL 20-100µL 50-200µL	2-20µL 10-100µL 20-200µL	2-20µL 20-100µL 50-200µL	2-20µL 20-100µL 50-200µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	100-1000µL 200-1000µL	1-5mL	1-10mL
2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	1-5mL	1-10mL
2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	2-20µL 10-100µL 20-200µL	100-1000µL	100-1000µL	100-1000µL	100-1000µL	1-5mL	1-10mL