

RPJ819Mu01 10µg

Recombinant Zinc Finger Homeobox Protein 1B (ZFHX1B)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Gln962~Glu1215

Tags: N-terminal His Tag

Subcellular Location: Nucleus

Purity: > 97%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 4.9

Predicted Molecular Mass: 33.2kDa

Accurate Molecular Mass: 38/33/28kDa as determined by SDS-PAGE reducing conditions.

Phenomenon explanation:

The possible reasons that the actual band size differs from the predicted are as follows:

1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affect the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

QGDLLDGAQ DYMSGLDDMT DSDSCLSRKK IKKTESGMYA
 CDLCDKTFQK SSSLLRHKYE HTGKRPHQCQ ICKKAFKHKH HLIHSRLHS
 GEKPYQCDKC GKRFSHSGSY SQHMNHRYSY CKREAEEEA AEREAREKGH
 LEPT E

[IDENTIFICATION]

XAGGGAGACTTGGTGGATGGGACAGGACTACATGTCAGGGCTAGATGACATGACAGACTCCGATTCCCTGCTGCTGAAAGAGATAAAGAGACAGAAAGTGGCATGTATGCGATGTGACTTATGTGCAAGACTTCCGAAAAGCAGTTCCCTTCCTCCGACATAAATAGGAACACACAGGAAAGAGACACACACAGTGTGGA
 Q G D L L D G A Q D Y M S G L D D M T D S D S C L S R K K I R K T E S G M Y A C D L C D R T F Q R S S L L R H K Y E H T G K R P H Q C Q I




Figure. Gene Sequencing (Extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.