

Horizontal Electrophoresis Apparatus and Power Packs

Hi-Eco Mini Horizontal Gel Box LA844


Features

- ◆ Robust acrylic frame mounted on non-slip rubber feet.
- ◆ A safety lid fits on one way only preventing access to live components.
- ◆ Fully shrouded connectors on safety lid.
- ◆ Strong thick removable gel casting tray resists warping.
- ◆ End gaskets allow gel pouring in the running tank or in an external casting unit.
- ◆ Long life pure platinum electrodes. Convenient carrying handles enable easy buffer disposal and increased mechanical strength.
- ◆ Operating temperature range 4-65°C.

Sr No	Name of Item	Specifications
1.	Principal Material	Acrylic
2.	Inner Tank Dimension L x W x H	16.2 x 9.9 x 2.9 cm
3.	No. of combs	2 Nos.
4.	No. of gel casting tray	1 no.
5.	Connecting cord	1set (black and 1Red)
6.	No. of platinum electrodes	1 anode 1 cathode
7.	Lid	1 No.
8.	Recommended power supply	LA842
9.	Gel Dimensions	10.0 X 7.5 cm
10.	Comb thickness	1.5mm polycarbonate
11.	Comb configurations	8 tooth, 8 tooth
12.	Voltage Limit	100 volts
13.	Max. Buffer volume	150 ml

Ordering Information

Hi-Eco Mini Horizontal Gel Box LA844



Mini Horizontal LA665 Electrophoresis System

Features

- ◆ Robust acrylic frame mounted on non-slip rubber feet.
- ◆ A safety lid fits on one way only preventing access to live components.
- ◆ Fully shrouded connectors on safety lid.
- ◆ Strong thick removable gel casting tray resists warping.
- ◆ End gaskets allow gel pouring in the running tank or in an external casting unit.
- ◆ Long life pure platinum electrodes.
- ◆ Convenient carrying handles enable easy buffer disposal and increased mechanical strength.
- ◆ Operating temperature range 4-65°C.

Sr No	Name of Item	Specifications
1.	Principal Material	Acrylic
2.	Inner Tank Dimension L x W x H	17.6 x 12 x 9.5cm
3.	No of combs	2 Nos.
4.	No. of gel casting tray	1 no.
5.	Connecting cord	1set (black and 1Red)
6.	No. of platinum electrodes	1 anode, 1 cathode
7.	Lid	1 No.
8.	Recommended power supply	LA690
9.	Gel Dimensions	7 x 8 cm (gasketed)
10.	Comb thickness	1.5mm polycarbonate
11.	Comb configurations	6 & 10 tooth
12.	Voltage Limit	150 volts
13.	Max. Buffer volume	400 ml

Ordering Information

Mini Horizontal LA665 Electrophoresis System

