



Horizontal Electrophoresis System LA666

Features


- ♦ Robust acrylic frame mounted on non-slip rubber feet.
- ♦ A safety lid fits on one way only preventing access to live components.
- ♦ Fully shrouded connectors on safety lid.
- ♦ Strong thick removable gel casting tray resists warping.
- ♦ End gaskets allow gel pouring in the running tank or in an external casting unit.
- ♦ Long life pure platinum electrodes.
- ♦ Convenient carrying handles enable easy buffer disposal and increased mechanical strength.
- ♦ Operating temperature range 4-65°C.

Sr No	Name of Item	Specifications
1.	Principal Material	Acrylic
2.	Inner Tank Dimension L x W x H	22 x 15 x 9.5cm
3.	No of combs	2 Nos.
4.	No. of gel casting tray	1 no.
5.	Connecting cord	1set (black and 1Red)
6.	No. of platinum electrodes	1 anode, 1 cathode
7.	Lid	1 No.
8.	Recommended power supply	LA690
9.	Gel Dimensions	9 x 11cm (gasketed)
10.	Comb thickness	1.5mm polycarbonate
11.	Comb configurations	10 & 14 tooth
12.	Voltage Limit	150 volts
13.	Max. Buffer volume`	600 ml

Ordering Information

**Horizontal
Electrophoresis System**

LA666



Horizontal Gel Box LA884

Features

- ♦ Robust acrylic frame mounted on non-slip rubber feet.
- ♦ A safety lid fits on one way only preventing access to live components.
- ♦ Fully shrouded connectors on safety lid.
- ♦ Strong thick removable gel casting tray resists warping.
- ♦ End gaskets allow gel pouring in the running tank or in an external casting unit.
- ♦ Long life pure platinum electrodes.
- ♦ Convenient carrying handles enable easy buffer disposal and increased mechanical strength.
- ♦ Operating temperature range 4-65°C.

Sr No	Name of Item	Specifications
1.	Principal Material	Acrylic
2.	Inner Tank Dimension L x W x H	24.5 x 18 x 9.5 cm
3.	No of combs	2 Nos.
4.	No. of gel casting tray	1 no.
5.	Connecting cord	1set (black and 1Red)
6.	No. of platinum electrodes	1 anode 1 cathode
7.	Lid	1 No.
8.	Recommended power supply	LA690
9.	Gel Dimensions	12x14cm (gasketed)
10.	Comb thickness	1.5mm polycarbonate
11.	Comb configurations	12 & 20 tooth
12.	Voltage Limit	150 volts
13.	Max. Buffer volume ®	800 ml

Ordering Information

Mini Horizontal Gel Box

LA884

