

## 100 well Multi channel compatible Horizontal Gel Box LA839


### Features

- ♦ Robust acrylic frame mounted on non-slip rubber feet.
- ♦ A safety lid fits on one way only preventing access to live components.
- ♦ Fully shrouded connectors on safety lid.
- ♦ Strong thick removable gel casting tray resists warping.
- ♦ End gaskets allow gel pouring in the running tank or in an external casting unit.
- ♦ Long life pure platinum electrodes.
- ♦ Convenient carrying handles enable easy buffer disposal and increased mechanical strength.
- ♦ Operating temperature range 4-65°C.

Sr No	Name of Item	Specifications
1.	Principal Material	Acrylic
2.	Inner Tank Dimension L x W x H	36 x 29.5 x 8cm
3.	No of combs	4 Nos.
4.	No. of gel casting tray	1 no.
5.	Connecting cord	1set (black and 1Red)
6.	No. of platinum electrodes	1anode 1 cathode
7.	Lid	1 No.
8.	Recommended power supply	LA690
9.	Gel Dimensions	25 x20 cm
10.	Comb thickness	1.5mm polycarbonate
11.	Comb configurations	25 tooth X 4 combs
12.	Voltage Limit	200 volts
13.	Max. Buffer volume	800 ml

### Ordering Information

**100 well Multi channel compatible Horizontal Gel Box LA839**



## Electrophoresis Power Supply LA842

### Features

- ♦ Two apparatus can be connected to the power supply with dual output at both side of the power supply.
- ♦ After switching ON the mains supply power, red led will glow.
- ♦ Select 50V 110 Volts DC voltage from selectable switch at front panel.

Sr No	Name of Item	Specifications
1.	Linear Power supply	50V/100 VDC Selectable, 200 mA
2.	Input	230 V AC
3.	Maximum operating temperature	45°C
4.	No. of output	2 nos
5.	Dimension in mm L x W x H	105 x 125 x 65
6.	Recommended for	LA844

### Ordering Information

**Electrophoresis Power Supply LA842  
(2 Terminal out put)**

