

Product Name	Code
Yeast and Fungal DNA Kits	
HiPurA™ Yeast Genomic DNA Purification Kit	MB552
HiPurA™ 96 Yeast Genomic DNA Purification Kit	MB555
HiPurA™ Yeast Plasmid DNA Purification Kit	MB550
HiPurA™ Fungal DNA Purification Kit	MB543
Environmental Sample DNA Kits	
HiPurA™ Soil DNA Purification Kit	MB542
HiPurA™ Water DNA Purification Kit	MB547
Mammalian DNA Kits	
HiPurA™ Mammalian Genomic DNA Purification Kit	MB506
HiPurA™ Sperm Genomic DNA Purification Kit	MB522
Forensic Sample DNA Kits	
HiPurA™ Forensic Sample Genomic DNA Purification Kit	MB524
HiPurA™ Bone DNA Purification Kit	MB525
HiPurA™ Buccal DNA Purification Kit	MB531
Agarose Gel DNA and PCR Product Purification Kits	
HiPurA™ Silica Kit for DNA isolation	MB503
HiPurA™ Quick Gel Purification Kit	MB539
HiPurA™ PCR Product Purification Kit	MB512
HiPurA™ 96 PCR Product Purification Kit	MB519
HiPurA™ PCR Product and Gel Purification Combo Kit	MB563
Food DNA Kits	
HiPurA™ Food DNA Purification Kit	MB562
HiPurA™ Food Pathogen DNA Purification Kit [Provided with culture media] New	MB568
HiPurA™ Food Pathogen (bacteria) DNA Purification Kit New	MB568A
Kits for Various Applications	
HiPurA™ Multi-Sample DNA Purification Kit	MB554
HiPurA™ Insect DNA Purification Kit	MB529
HiPurA™ Stool DNA Purification Kit	MB544
RDP Trio™ Reagent	MB566
HiPurA™ Paraffin-Embedded Tissue DNA Purification Kit	MB530
HiPurA™ Spice DNA Purification Kit	MB567
HiPurA™ White Spot Syndrome Virus (WSSV) DNA Purification Kit New	MB572
RNA Isolation Kits	
RNA-XPress™ Reagent	MB601
HiPurA™ Total RNA Miniprep Purification Kit	MB602
HiPurA™ 96 Total RNA Purification Kit	MB612
HiPurA™ Plant and Fungal RNA Miniprep Purification Kit	MB603
HiPurA™ Yeast RNA Purification Kit	MB611
HiPurA™ Bacterial RNA Purification Kit	MB613
HiPurA™ Viral RNA Purification Kit* New	MB615



Product Name	Code
HiPurA™ All Blood RNA Purification Kit	MB617

Molecular Biology Buffers/Reagents for DNA/

RNA Isolation

RNase A Solution (20 mg/ml)	DS0003
Proteinase K Solution (20 mg/ml)	DS0013
Lysis Solution (C1)	DS0010
Pre Wash Solution Concentrate (PW)	DS0011
Wash Solution Concentrate (WS)	DS0012
Lysis Solution (AL)	DS0015
Resuspension Solution (HP1)	DS0020
Lysis Solution (HP2)	DS0021
RNA Lysis Solution (HRL) New	DS0037
Neutralization Solution (HN3)	DS0022
Wash Solution Concentrate (HPE)	DS0024
PCR Binding Solution (SPB)	DS0027
Wash Solution (HPB)	DS0032
Elution Buffer (ET) (10 mM Tris-Cl, pH 8.5)	DS0040
10X RBC Lysis Buffer	R075
CTAB Extraction Solution	ML066
DNase I Solution (1mg/ml)	ML068

InstaDNA™ Extraction

InstaDNA™ kit	MBT109
---------------	--------

InstaDNA™ Cards

InstaDNA™ card	MBT126
InstaDNA™ Twin Card New	MBT137
InstaDNA™ Quadricard	MBT121

Ladders/Markers

DNA Ladders

50bp DNA ladder	MBT084
100bp DNA ladder	MBT049
100bp DNA ladder M	MBT130
1kb DNA ladder	MBT051
High Range DNA Ladder	MBT090

Protein Ladder

Prestained Protein Ladder (10 - 245 kDa)	MBT092
--	--------