

## Nucleic Acid Amplification

### PCR KITS FOR DIAGNOSTIC APPLICATIONS

#### Kits For Clinical Diagnostics

##### Real-Time PCR Kits

(Provided With Positive Control)

###### Probe Based Kits

Salmonella Detection Kit (Real-Time Probe Based PCR) <b>New</b>	MBPCR099
Generic Dengue Detection Kit One -Step (Real-Time Probe Based PCR) <b>New</b>	MBPCR101
HPV Detection Kit (Real-Time Probe Based PCR) <b>New</b>	MBPCR105
Mycobacterium tuberculosis Detection Kit (Real-Time Probe Based PCR) <b>New</b>	MBPCR108
Malaria Detection Kit (Real-Time Probe Based PCR) <b>New</b>	MBPCR111
Chikungunya Detection Kit (Real-Time Probe Based PCR) <b>New</b>	MBPCR112

###### SYBr Based Kits

Mycoplasma Detection Kit	MBPCR015
Mycobacterium tuberculosis Detection Kit	MBPCR017
Multi-Drug Resistant Mycobacterium tuberculosis Detection Kit (Uniplex) <b>New</b>	MBPCR073
E. coli O157:H7 Detection Kit	MBPCR023
Salmonella Detection Kit	MBPCR024
Listeria monocytogenes Detection Kit	MBPCR025
Campylobacter jejuni Detection Kit	MBPCR026
Vibrio cholerae Detection Kit	MBPCR027
Legionella pneumophila Detection Kit	MBPCR028
Chronobacter sakazakii Detection Kit	MBPCR029
Klebsiella pneumoniae Detection Kit	MBPCR057
Pseudomonas aeruginosa Detection Kit	MBPCR058
Generic E.coli Detection Kit	MBPCR059
Candida albicans Detection Kit	MBPCR060
Bacillus subtilis Detection Kit	MBPCR061
Shigella Spp. Detection Kit	MBPCR062
Staphylococcus aureus Detection Kit	MBPCR064
Acinetobacter baumannii Detection kit	MBPCR065
Malaria Detection Kit (Real-Time) <b>New</b>	MBPCR092
White Spot Syndrome Virus (WSSV) Detection Kit <b>New</b>	MBPCR093
Generic Dengue Detection Kit (One-Step Real Time SYBr Based PCR) <b>New</b>	MBPCR103
Chikungunya Detection Kit (One Step Real time PCR)	MBPCR083

##### Semi-Quantitative PCR Kits

(Provided With Positive & Internal Controls)

Campylobacter jejuni Detection Kit	MBPCR001
E. coli O157.H7 Detection Kit	MBPCR002
Vibrio cholerae Detection Kit	MBPCR003
Salmonella Detection Kit	MBPCR004
Listeria monocytogenes Detection Kit	MBPCR005
Legionella pneumophilla Detection Kit	MBPCR006
White Spot Syndrome Virus (WSSV) Detection Kit <b>New</b>	MBPCR013
Methicillin Resistant Staphylococcus aureus (MRSA) Detection Kit (Uniplex)	MBPCR008
Methicillin Resistant Staphylococcus aureus (MRSA) Detection Kit (Multiplex)	MBPCR020
Staphylococcus Cassette Chromosome (SCC) mec typing MRSA Detection Kit (Multiplex)	MBPCR019
Mycobacterium tuberculosis Detection Kit	MBPCR009
Multi-Drug Resistant Mycobacterium tuberculosis Detection Kit (Multiplex)	MBPCR022
Malaria Detection Kit (Multiplex)	MBPCR016
Chronobacter sakazakii Detection Kit	MBPCR030
Acinetobacter baumannii Detection kit	MBPCR038
Diarrhoeogenic E.coli Detection Kit (Multiplex)	MBPCR039
Klebsiella pneumoniae Detection Kit	MBPCR044
Pseudomonas aeruginosa Detection Kit	MBPCR050
Generic E.coli Detection Kit	MBPCR051
Candida albicans Detection Kit	MBPCR052
Bacillus subtilis Detection Kit	MBPCR053
Shigella spp. Detection Kit	MBPCR054
Staphylococcus aureus Detection Kit	MBPCR056
Clostridium sporogenes Detection Kit <b>New</b>	MBPCR074
Beta-Thalassemia Detection Kit (Multiplex)	MBPCR014
Beta-Thalassemia Detection Kit (Uniplex) <b>New</b>	MBPCR014A
Generic Dengue Detection Kit (One-Step Semi-Q PCR Based) <b>New</b>	MBPCR102
Chikungunya Detection Kit <b>New</b> (One Step Semi quantitative PCR)	MBPCR084
Chikungunya Detection Kit (Two Step Semi quantitative PCR)	MBPCR085