



Product Name	Code
Palm PCR Kits	
Palm E.coli O157:H7 Detection Kit	MBPCR032
Palm Salmonella Detection Kit	MBPCR033
Palm Mycobacterium tuberculosis Detection Kit	MBPCR034
Palm Malaria Detection Kit	MBPCR035
Palm Listeria monocytogenes Detection Kit	MBPCR036
Palm Multi Drug Resistant Mycobacterium tuberculosis (MDR-TB) Detection Kit	MBPCR037
Palm Bacillus subtilis Detection Kit	MBPCR040
Palm Candida albicans Detection Kit	MBPCR041
Palm Pseudomonas aeruginosa Detection Kit	MBPCR042
Palm Staphylococcus aureus Detection Kit	MBPCR043
Palm Generic E.coli Detection Kit	MBPCR045
Palm Diarrhoeagenic E.coli Detection Kit	MBPCR046
Palm Shigella Spp. Detection Kit	MBPCR047
Palm GMO (Genetically Modified Organisms) Detection Kit	MBPCR048
Palm Klebsiella pneumoniae Detection Kit	MBPCR049
Palm Acinetobacter baumannii Detection kit	MBPCR066
Palm Clostridium sporogenes Detection Kit New	MBPCR067
Palm Salmonella enterica subsp. enterica Detection Kit New	MBPCR068
Palm Aspergillus brasiliensis Detection Kit New	MBPCR069
Palm Generic Dengue Detection Kit (One Step RT-PCR) New	MBPCR082
Palm White Spot Syndrome Virus (WSSV) Detection Kit New	MBPCR100
Palm Chikungunya Detection Kit New	MBPCR086

Kits For Food Diagnostics

Extraction Kits	
HiPurA™ Food DNA Purification Kit	MB562
HiPurA™ Food Pathogen DNA Purification Kit (Salmonella) (Provided with culture media) New	MB568
HiPurA™ Food Pathogen (Bacteria) DNA Purification Kit New	MB568A

Food PCR Kits	
Real Time PCR Kits (With Positive Control) Probe Based Food PCR Kits	
Salmonella Food Detection Kit (Real-time PCR Kit Probe Based) New	MBPCR099

Product Name	Code
SYBr Based Food PCR Kits	
Generic E. coli Food Detection Kit (Real-time SYBr Based PCR Kit) New	MBPCR059A
Salmonella Food Detection Kit (Real-time SYBr Based PCR Kit) New	MBPCR024A
Listeria monocytogenes Food Detection Kit (Real-time SYBr Based PCR Kit) New	MBPCR025

Semi-Quantitative PCR Kits (Provided with Positive & Internal Controls)	
Salmonella Food Detection Kit (Semi-Quantitative PCR Kit) New	MBPCR004A

Food Extraction-PCR Combo Kits	
E.coli O157:H7 and Salmonella spp. Food Detection Kit	MBT112
E.coli O157:H7 and Salmonella spp. Food Detection Kit (Provided with culture media)	MBT113

Food Heat lysis-PCR Combo Kits	
E.coli O157:H7 Food Detection Kit	MBT114
E.coli O157:H7 Food Detection Kit (Provided with culture media)	MBT115
Salmonella Food Detection Kit	MBT116
Salmonella Food Detection Kit (Provided with culture media)	MBT117

Meat Adulteration/Identification Kits	
Pork Detection Kit (Real-time SYBr based PCR Kit) New	MBPCR098

Kits For GMO Detection	
GMO (Genetically Modified Organism) Detection Kit (Semi-Q PCR Based)	MBPCR055
GMO (Genetically Modified Organism) Detection Kit (Real-time PCR Based)	MBPCR063