



Product Name

Code

Product Name

Code

POSITIVE CONTROLS FOR PCR

Pseudomonas aeruginosa Genomic DNA	MBT097	Listeria monocytogenes Genomic DNA	MBT104
Klebsiella pneumoniae Genomic DNA	MBT098	Legionella pneumophila Genomic DNA	MBT105
Shigella flexneri Genomic DNA	MBT099	Campylobacter jejuni Genomic DNA	MBT106
Streptococcus pneumoniae Genomic DNA	MBT100	Mycobacterium tuberculosis Genomic DNA	MBT107
Vibrio cholerae Genomic DNA	MBT101	Acinetobacter baumannii Genomic DNA	MBT122
E. coli O157:H7 Genomic DNA	MBT102	Campylobacter coli Genomic DNA New	MBT135
Salmonella Genomic DNA	MBT103		

Cloning Reagents

Restriction Enzymes

EcoR I	MBRE001
BamH I	MBRE002
Hind III	MBRE003
Taq I	MBRE004
Not I	MBRE005
Hinf I	MBRE006
Hae III	MBRE007
Alu I	MBRE008
EcoR V	MBRE009
Kpn I	MBRE010
Pst I	MBRE011
Sac I	MBRE012
Sca I	MBRE013
Sph I	MBRE014
Stu I	MBRE015
Xba I	MBRE016
Nco I	MBRE017
Bgl II	MBRE018
Hinc II	MBRE019
Mlu I	MBRE020
Sma I	MBRE021
Pvu II	MBRE022

Modifying Enzymes

T4 DNA Ligase (200U/ml)	MBT091
-------------------------	--------

Competent Cells

HiPurA™ Competent Cells	MBT111
-------------------------	--------

Transformation Kit

Bacterial Transformation kit	MBT127
------------------------------	--------

Protein Purification and Crystallography

Protein Purification

IgG Purification Kit (Protein A Based)	
IgG Purification Kit (Protein A Based, Gravity Flow)	MBPP001
IgG Purification Kit (Protein A Based, Spin Column)	MBPP001SP
IgG Purification Kit (Protein G Based)	
IgG Purification Kit (Protein G Based, Gravity Flow)	MBPP002
IgG Purification Kit (Protein G Based, Spin Column)	MBPP002SP

IgG Purification (Protein L Based)

IgG Purification Kit (Protein L Based, Gravity Flow)	MBPP004
IgG Purification Kit (Protein L Based, Spin Column)	MBPP004SP

GST-Tagged Bacterial Protein Purification Kit

GST-Tagged Bacterial Protein Purification Kit (Gravity Flow)	MBPP003
GST-Tagged Bacterial Protein Purification Kit (Spin Column)	MBPP003SP

HIS Tagged Bacterial Protein Purification Kit

HIS-Tagged Bacterial Protein Purification Kit (Gravity Flow)	MBPP005
HIS-Tagged Bacterial Protein Purification Kit (Spin Column)	MBPP005SP

Protein Crystallography

Protein Crystallization Kit	MBPC001
Protein Standards	MBT129