

Product Name	Code
HiPer® Bacterial Conjugation Teaching Kit	HTM004
HiPer® Transduction Teaching Kit	HTM005
HiPer® Ezy MIC™ Teaching Kit	HTM006
HiPer® Protein Assay Teaching Kits	
HiPer® SDS-PAGE Teaching Kit	HTP001
HiPer® Protein Silver Staining Teaching Kit	HTP002
HiPer® Native PAGE Teaching Kit	HTP003
HiPer® Immunology Teaching Kits	
HiPer® Latex Agglutination Teaching Kit	HTI001
HiPer® Ouchterlony Double Diffusion Teaching Kit (For Antigen-Antibody Pattern)	HTI002
HiPer® Ouchterlony Double Diffusion Teaching Kit (For Antibody Titration)	HTI003
HiPer® Radial Immunodiffusion Teaching Kit	HTI004
HiPer® Immunoelectrophoresis Teaching Kit	HTI005
HiPer® Rocket Immunoelectrophoresis Teaching Kit	HTI006
HiPer® Counter Current Immunoelectrophoresis Teaching Kit	HTI007
HiPer® Blood Grouping Teaching Kit	HTI008
HiPer® Western Blotting Teaching Kit	HTI009
HiPer® Quantitative Precipitin Assay Teaching Kit	HTI011
HiPer® Antigen Capture ELISA Teaching Kit	HTI012
HiPer® Antibody Capture ELISA Teaching Kit	HTI013
HiPer® Sandwich ELISA Teaching Kit	HTI014
HiPer® Dot ELISA Teaching Kit	HTI015
HiPer® Immunoprecipitation Teaching Kit	HTI016
HiPer® Widal Test Teaching Kit (Slide Test)	HTI017
HiPer® Widal Test Teaching Kit (Tube Test)	HTI018
HiPer® RA Test Teaching Kit	HTI019

Product Name	Code
HiPer® RPR Test Teaching kit	HTI020
HiPer® VDRL Test Teaching Kit	HTI021
HiPer® IgG Purification Teaching Kit	HTI022
HiPer® Biochemistry Teaching Kits	
HiPer® Yeast Cell Immobilization Teaching Kit	HTBC001
HiPer® Carbohydrates Estimation Teaching Kit (Qualitative)	HTBC002
HiPer® Carbohydrates Estimation Teaching Kit (Quantitative)	HTBC003
HiPer® Protein Estimation Teaching Kit (Qualitative)	HTBC004
HiPer® Protein Estimation Teaching Kit (Quantitative)	HTBC005
HiPer® DNA Estimation Teaching Kit	HTBC006
HiPer® RNA Estimation Teaching Kit	HTBC007
HiPer® Chromatography Teaching Kits	
HiPer® Ion Exchange Chromatography Teaching Kit	HTC001
HiPer® Affinity Chromatography Teaching Kit	HTC002
HiPer® Thin Layer Chromatography Teaching Kit	HTC003
HiPer® Gel Filtration Chromatography Teaching Kit	HTC004

Latex Agglutination Test Kits

HiClostridium difficile Latex Test Kit	LK01
HiSalmonella™ Latex Test Kit	LK02
HiStaph Latex Test Kit	LK03
HiLegionella Latex Test Kit	LK04
HiStrep™ Latex Test Kit	LK06
HiListeria™ Latex Test Kit	LK07
HiRotavirus Latex Test Kit	LK08
HiE.coli™ Latex Test Kit	LK13

