



UNBOX RNA





Comparative PCR data of Dengue viral RNA purified from human tissue sample using HiMedia's HiPurA™ Viral RNA Purification Kit

Comparative data representing PCR amplification of cells RNA using HiMedia's HiScript Two Step RT-PCR kit (MBT138)




Comparative data of Dengue detection using HiMedia's HiScript One Step RT-PCR Kit (MBT095)

Comparative data of RNase Inhibitor (ML162)



Comparative data of RNA, extracted from rice leaves, preserved in RNA stabilizer (ML161) at room temperature for one week

RNA Stabilizer ML161



Chikungunya Detection using Hi-Quanti One Step RT-PCR Kit (MBT128)

