

Lab Aids

Universal tips

a) General Tips :

i) CG308 - 0.5 - 10 μ l ii) CG309 - 10 μ l iii) CG311 - 5000 μ l iv) LA888 - 10000 μ l

b) Micropipette Tips :

i) LA972 - 200 μ l ii) LA973 - 1000 μ l iii) LA975 - 1000 μ l

- ▶ Nuclease free
- ▶ Pyrogen free
- ▶ Autoclavable

c) Barrier Tips :




i) LA749 - 10 μ l ii) LA750 - 20 μ l iii) LA751 - 200 μ l iv) LA859 - 1000 μ l,

- ▶ Ultra low retention
- ▶ Pre sterile
- ▶ Nuclease & Endotoxin free

d) Non Barrier Tips :

i) LA788 - 10 μ l ii) LA789 - 200 μ l iii) LA790 - 1000 μ l

- ▶ Ultra low retention
- ▶ Pre sterile
- ▶ Nuclease & Endotoxin free



μ Pet™ Autoclavable Micropipette

a) Single channel Pipette : variable & fixed volume

i) LA955 - 0.5 - 10 μ l ii) LA959 - 20 - 200 μ l
 iii) LA961 - 500 - 5000 μ l iv) LA963 - 1000 μ l (fixed volume)

b) Multi channel Pipette : 8 channel

i) LA964 - 0.5 - 10 μ l ii) LA966 - 20 - 200 μ l iii) LA967 - 30 - 300 μ l

