

Accessories for Teaching Kits

Sr. no.	Plasticware/Glassware/Chemicals	Product code	Sr. no.	Plasticware/Glassware/Chemicals	Product code
1	Prima-DUO Thermal Cycler	LA948	16	15ml Centrifuge Tubes	PW1228
2	Prima-96 Thermal Cycler	LA949	17	Autoclavable L- Spreader	PW1085
3	Horizontal Electrophoresis System	LA666/LA844	18	Microscopic Slides	CG001
4	Electrophoresis Power supply	LA690/LA842	19	Microscope cover glasses (20X20 mm)	CG120
5	Micropipettes(Varivol)		20	HiGas Burner (With LPG connector for gas)	LA729
a	Micropipette-10 (0.5 to 10µl)	LA611	21	Inoculation Loops - 1µl	CG171
b	Micropipette-50 (5 to 50µl)	LA612	22	Test tube racks(Z - Shape)	CG260, CG261
c	Micropipette-200 (50 to 200µl)	LA613	23	Autoclavable Petri plates	PW008
d	Micropipette-1000 (200 to 1000µl)	LA614	24	Autoclavable Test Tube, Rimless	PW1162
e	Micropipette-100 (10 to 100µl)	LA615	25	Stainless steel forceps, Blunt and Pointed	LA706 - LA111
f	Micropipette-1000 (100 to 1000µl)	LA616	26	Spatula - 128	LA007
g	Micropipette-20 (2 to 20µl)	LA617	27	Ethidium Bromide (10mg/ml)	MB074
h	Micropipette-200 (20 to 200µl)	LA618	28	Molecular Biology Grade Water	MI024
6	Insta DNA™ card punch 3mm	LA1000	29	Sulphuric Acid, Ultra Grade	RM6244
7	Latex Examination Gloves (sterile & nonsterile), available in three size	LA891/LA893	30	Hydrochloric Acid, 35% Pure, A.R.	AS004
8	a Empty Tip box (For 200µl tips)	LA852	31	Methanol	MB113
b Empty Tip box (For 1000µl tips)	LA977		32	2- Mercaptoethanol	MB041
9	Hi-Pette Autopipette Stand	LA632	33	Chloroform	MB109
10	a Micropipettes Tips- 10µl	CG308	34	1-Octanol	GRM2933
b Micropipettes Tips- 200µl	CG304		35	Diluent for DNA Extraction	MB228
c Micropipettes Tips- 1000µl	CG305		36	Nitric Acid about 70% pure, L.R.	AS008
11	Barrier Tips- 10µl	LA749	37	DEPC (Diethyl Pyrocarbonate)	MB076
12	Barrier Tips- 20µl	LA750	38	Acetic Acid	MB052
13	Barrier Tips- 200µl	LA751	39	Hi-UV™ Transilluminator New	LA1047
14	Barrier Tips- 1000µl	LA859	40	TabSpin™ Microcentrifuge New	LA1001
15	50 ml Centrifuge Tubes	PW1230	41	InstaQ™ 48M4	LA023
			42	InstaQ 96™	LA1012



LA844



LA842



LA1001


TabSpin

LA1047



LA948

Prima-DUO™ Thermal Cycler



LA666

Mini Horizontal Gel Box



LA1012

Insta Q96



HiUV

