


Kits for Food Diagnostics (With Positive Control)

Probe Based Food PCR Kits

MBPCR099 Salmonella Detection Kit
(Real-time PCR Kit Probe Based)

Sample	C _t value
Positive control	14.89
<i>Salmonella</i> spiked in milk	15.32
Negative control	N/A

Real Time Probe Based PCR
for *Salmonella* Detection in Spiked Milk
(Code : MBPCR099)




SYBr Based Food PCR Kits

MBPCR025A *Listeria monocytogenes* Food Detection
Kit (Real-time SYBr Based PCR Kit)

Sample	C _t value
Positive control	9.96
<i>L. monocytogenes</i> spiked in milk	10.85, 10.81
Negative control	N/A

Real Time SYBr Based PCR Data
for *L. monocytogenes* Detection in Spiked Cheese
(Code : MBPCR025A)



Meat Adulteration / Identification Kits

MBPCR136 Pork Detection Kit
(Real-Time Probe based PCR Kit)

Sample	C _t value
Pork DNA	25.81
Internal Control	17.92

Pork Detection
(Sensitivity : 0.01%)

