


MBPCR138 Buffalo Detection Kit
(Real-Time Probe based PCR)


Data for Buffalo Detection Kit
(MBPCR138) (Sensitivity : 0.5%)



Sample	C _t value
Buffalo milk in cow milk	19.02
Positive control	15.15
Negative control	N/A

MBPCR139 Cow Detection Kit
(Real-Time Probe based PCR)

Data for Cow Detection Kit
(MBPCR139)




Sample	C _t value
Cow milk in Buffalo milk	21.2
Positive control	16.73
Negative control	N/A

Kits for GMO Detection

MBPCR063 GMO (Genetically Modified Organism)
Detection Kit (Real-time PCR Based)

Gene	C _t value
35S	15.43, 15.98
NOS	20.19, 19.89
rbcl	13.9, 13.4
Negative Control	N/A

Data representing GMO detection
in cotton plant



Kits for Veterinary Diagnostics (With Positive Control)

MBPCR129 Animal Hemorrhagic septicemia (HS)
Detection Kit (Real-Time SYBr Based PCR Kit)

Sample	C _t value
Negative control	N/A
Positive control	09.45

Data for HS Detection (MBPCR129)

