

Ordering Information

Product Description	Kit Components	Product Code
47mm Filtration Assembly		
SF198 series (Autoclavable, Stainless Steel) includes <ul style="list-style-type: none"> • Filter holder : 47mm, 250ml • #8 size rubber stop cork • duck clamp 	1000ml receiver flask	SF198
	1000ml receiver flask w/ Vacuum pump	SF198A
	1000ml receiver flask w/ Vacuum trap assembly	SF198B
	1000ml receiver flask w/ Vacuum trap assembly & vacuum pump	SF198C
SF199 series (Autoclavable, Glass) includes <ul style="list-style-type: none"> • Filter holder : 47mm, 250ml • #8 size rubber stop cork • duck clamp 	1000ml receiver flask	SF199
	1000ml receiver flask w/ Vacuum pump	SF199A
	1000ml receiver flask w/ Vacuum trap assembly	SF199B
	1000ml receiver flask w/ Vacuum trap assembly & vacuum pump	SF199C
SF200 series (Autoclavable, Polypropylene) includes <ul style="list-style-type: none"> • Filter holder : 47mm, 250ml • #8 size rubber stop cork • duck clamp 	1000ml receiver flask	SF200
	1000ml receiver flask w/ Vacuum pump	SF200A
	1000ml receiver flask w/ Vacuum trap assembly	SF200B
	1000ml receiver flask w/ Vacuum trap assembly & vacuum pump	SF200C
142mm Filtration Assembly		
SF202 series (Autoclavable, Stainless Steel) includes <ul style="list-style-type: none"> • Filter holder : 142mm 	1000ml receiver flask	SF202
	1000ml receiver flask w/ Vacuum pump	SF202A
	1000ml receiver flask w/ Vacuum trap assembly	SF202B
	1000ml receiver flask w/ Vacuum trap assembly & vacuum pump	SF202C
S.S. Sterility Test Manifold Assembly		
SF85 series (Autoclavable, Stainless Steel) includes <ul style="list-style-type: none"> • 3 Stations 	1000ml receiver flask	SF85
	1000ml receiver flask w/ Vacuum pump	SF85A
	1000ml receiver flask w/ Vacuum trap assembly	SF85B
	1000ml receiver flask w/ Vacuum trap assembly & vacuum pump	SF85C
SF86 series (Autoclavable, Stainless Steel) includes <ul style="list-style-type: none"> • 6 Stations 	1000ml receiver flask	SF86
	1000ml receiver flask w/ Vacuum pump	SF86A
	1000ml receiver flask w/ Vacuum trap assembly	SF86B
	1000ml receiver flask w/ Vacuum trap assembly & vacuum pump	SF86C

Note : All assemblies available with 2000ml and 5000ml receiver flask.



For Life is Precious

HiMediaLaboratories™
www.himedialabs.com

CORPORATE OFFICE

A-516, Swastik Disha Business Park, Via Vadhani Indl Est, LBS Marg, Mumbai - 400 086, India.

Tel : +91-22-6147 1919 / 2500 3747 | Fax : +91-22-6147 1920 / 2500 5764

Email : atc@himedialabs.com

