


QC Microbiology



QC-Microbiology			
		0.22 μ	0.45 μ
Sterile Membranes	Cellulose nitrate (47mm)	-	SF95H
	Gridded Cellulose nitrate (47mm)	SF97B	SF97D
	MCE (47mm)	SF99A	SF99B
Sampling Manifold System	3 Stations	S.S. Sterility Test Manifold System (SF85) with 3 Stations & without holders	
	6 Stations	S.S. Sterility Test Manifold System (SF86) with 6 Stations & without holders	
Filtration Assembly	142 mm	Stainless Steel Filtration Assembly (SF202A/B/C) Contents: 142mm Filter holder Receiver Flasks (1000mL/2000mL/5000mL)	
Filter Holders	13mm	S.S. Sterility test Filter Holder for 15ml (SF80)	
	25 mm	S.S. Sterility test Filter Holder for 150ml (SF81)	
	47 mm	S.S. Sterility test Filter Holder for 300ml (SF82)	
	47 mm	S.S. Sterility test Filter Holder for 500ml (SF83)	
	47 mm	S.S. Sterility test Filter Holder for 1000ml (SF84)	
Filtration Funnel	With Membrane	Filtration Funnel – SF191 with stop cork	
	Without Membrane	Filtration Funnel – SF190 with 0.45μm edge hydrophobic Cellulose Nitrate membrane & without stop cork	
Vacuum Pump	-	Vacuum Pump with Hose Tube – SF201	