


QC Analytical



QC Analytical Labs (Non-Sterile Application)

| | | <i>with Prefilter</i> | | <i>without Prefilter</i> | | |
|----------------------------|------------------------------------|---|---------------|--|---|--|
| | | 0.22 μ | 0.45 μ | 0.22 μ | 0.45 μ | |
| Syringe Filters | Nylon-66 (25mm) | SF128 | SF129 | SF126 | - | |
| | (30mm) | - | - | SF152 | SF153 | |
| | PVDF (25mm) | SF130 | SF131 | - | - | |
| | (30mm) | - | - | SF154 | SF155 | |
| | PTFE (25mm) | SF132 | SF133 | - | - | |
| | (30mm) | - | - | SF156 | SF157 | |
| Membranes | | Nylon-66 | | PVDF | | |
| | | 47mm | 142mm | 47mm | 142mm | |
| | 0.22 μ | SF103A | SF105A | SF182A | SF183A | |
| | 0.45 μ | SF103B | SF105B | SF182B | SF183B | |
| | 0.8 μ | SF103D | SF105C | - | - | |
| | 3.0 μ | SF103H | - | SF182C | SF183C | |
| | 5.0 μ | SF103J | - | - | - | |
| | | HPLC Grade Nylon-66 Membranes - 47mm | | | | |
| | | 0.22 μ | 0.45 μ | 0.8 μ | All assembly other than those mentioned here are available on request | |
| | | SF103E | SF103F | SF103G | | |
| Filtration Assembly | 47mm | Solvent Filtration Assembly (Stainless steel – SF198, Glass – SF199, Polypropylene – SF200) | | | | |
| | 142 mm | Solvent Filtration Assembly (142mm/293mm Filter holder) | | | | |
| Filter Holders | | without integrity tester | | with integrity tester | | |
| | 90mm | S.S. Inline Process Filter Holder (SF87) | | S.S. Inline Process Filter Holder (SF88) | | |
| | 142mm | S.S. Inline Process Filter Holder (SF89) | | S.S. Inline Process Filter Holder (SF90) | | |
| | 293mm | S.S. Inline Process Filter Holder (SF91) | | S.S. Inline Process Filter Holder (SF92) | | |
| Vacuum Pump | Vacuum Pump with Hose Tube – SF201 | | | | | |