

# HiAnaerobic™ WORK STATION

Refer Page  
**684**

Smallest Bench Top  
Work Station for  
Growing Anaerobes




(LA1022)

## PERFECT CHOICE OF EASY TO USE ANAEROBIC AND MICROAEROPHILIC INCUBATED WORK STATIONS

- Ergonomically designed
- Ezee sleeve system
- Easy loading and unloading of work into the chamber
- Economical dual gas operation
- Exceptionally low gas consumption
- 6 mm gauge acrylic with 15mm laminated end panels
- 800 mm (W) x 660 mm (D) x 650 mm (H)
- Capacity to incubate upto 270 mm x 90 mm plates
- Operating temperature range ambient + 5°C to 45°C



EZEE SLEEVE SYSTEM



ANAEROBIC INCUBATION



ANAEROBIC INDICATOR



## LAB AIDS AND INSTRUMENTATION

Anaerobic System Mark II, Anaerobic System Mark III, Anaerobic System Mark V, Anaerobic System Mark VI, Anaero Indicator Tablet, Anaerogas Pack 3.5L/1.5L, Air Purifiers, Air Sampler Systems, Air Petri Sampling Systems, HiBio-ID™ Reader, Antimicrobial System & Antibiotic Zonescale, Colony Counter, Automatic Loop Sterilizer, Flat Bottom pH Electrode, Plate Master for Speed Plating, Needle & Syringe destroyer, Auto Smear Systems, Automatic Staining Systems, Plate Inspection Port.

[www.himedialabs.com](http://www.himedialabs.com)