

## Durham Tubes

**GW163**

Durham Tubes are very small size test tubes made from Neutral Glass. The Durham Tube is inverted tube inside the fermentation tube which actually captures the gas produced by microorganisms.

**Application :** Chemical laboratory, Biological laboratory , Microbiology & Diagnostic sectors , Industrial laboratory and various other laboratories.

Product Name	Product Code	Description	Size (mm)
Durham Tubes	GW163	These tube are made of neutral Glass. and they are autoclavable.	Length = 25-27 Diameter = 6-7

### Product features :

- Neutral Glass.
- Glass is very clear and transparent.
- The edges are properly polished and hence smooth.
- The bottoms of these Durham tubes are rounded.
- Pack Size : 1 X 100 Nos.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.