

MBT075

Hi-PCR Kit

Product Name	Product Code	Kit Packing**
Hi-PCR Kit	MBT075-10R	10 reactions
Hi-PCR Kit	MBT075-25R	25 reactions

**** The product is supplied with a vial of Molecular Biology Grade Water.**

Description:

Hi-PCR Kit is provided at affordable price and is designed to meet various needs while providing convenience in many ways. Hi-PCR kit contains all necessary components for PCR. The user needs to add the desired template and primers of interest. This kit is beneficial when sharing the reagents with other users.

Hi-PCR kit is provided with:

Components	Recommended volume to be added for 50 µl reaction	Final concentration in the reaction mixture	Reagents provided for 10R (reactions)	Reagents provided for 25R (reactions)
10X Assay Buffer	5 µl	1X	75 µl	150 µl
10 mM dNTP mix	5 µl	1 mM	75 µl	150 µl
50 mM MgCl ₂	2.5 µl	2.5 mM	75 µl	150 µl
Taq DNA Polymerase (1 U/µl)	1 µl	1 Units	50 µl	100 µl
Nuclease free water	Upto 50 µl	-	1 ml	2 ml
6X Loading Dye	-	-	100 µl	200 µl

Storage conditions: Store at -20°C. When stored under recommended conditions, the kit components are stable for 6 months.

Quality control:

Detected free of RNases, endonuclease and exonuclease activities.

Technical Assistance

At HiMedia, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.

PIMBT075_0/0914

MBT075-02

Disclaimer :

User must ensure of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. Himedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

CE A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai - 400 086, India Customer Care No.: 022-6147 1919 Email: techhelp@himedialabs.com



Registered Office :

23, Vadhani Industrial Estate, LBS Marg, Mumbai - 400 086, India.
Tel. : (022) 4017 9797 / 2500 1607
Fax : (022) 2500 2286

Commercial Office

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai - 400 086, India
Tel: 00-91-22-6147 1919
Fax: 6147 1920, 2500 5764
Email : info@himedialabs.com
Web : www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein