



**Product Information**

**Product Code : SF86-1NO**  
**S.S. Sterility Test Manifold System**  
**With 6 stations**  
**Without holders**

**Product description :**

S.S. Sterility Test Manifold System with 6 stations is an integrated filtration system having a prime application in filtering multiple samples off one vacuum source. It is without filter holders. It is perfectly supported on stainless steel manifold with rugged side supports which keeps the system rigid and firm even while fully loaded.

**Applications:**

- Sterility testing of antibiotics and pharmaceuticals.
- Vacuum filtration of biomass, waste water and particulate contamination analysis.
- Microbiological analysis of retained particles.
- Yeast and mold detection.
- Photoresists
- Vacuum filtration of solvents and aqueous solutions.

**Specifications:**

1. Construction material : Stainless steel 316, Stainless Steel 304 support plates, Teflon spool
2. Suitable filters : Sterility test filter holder - 47mm
3. Pressure : Vacuum only
4. Dimensions : 790mm X 150mm X 165mm
5. Weight : Approx 5Kg

**Packaging:**

Product	Code	Packing
S.S. Sterility Test Manifold system With 6 stations Without holders	SF86-1NO	1No