

Product	Code	Packing
5-Bromo-2'-deoxyuridine (BrdU) 5-BrdU Cell Culture Tested $C_9H_{11}BrN_2O_5$ MW : 307.10 CAS : 59-14-3 Assay : $\geq 96.5\%$ Store at 2 - 8°C	TC372-100MG TC372-250MG TC372-500MG	100mg 250mg 500mg
Cytidine Cell Culture Tested $C_9H_{13}N_3O_5$ MW : 243.22 CAS : 65-46-3 Assay : $\geq 99\%$ Store below 30°C	TC110-1G TC110-5G TC110-50G	1gm 5gm 50gm
2'-Deoxyadenosine monohydrate Cell Culture Tested $C_{10}H_{13}N_5O_3H_2O$ MW : 269.26 CAS : 16373-93-6 Assay : $\geq 99\%$ Store at 2 - 8°C	TC111-1G TC111-5G	1gm 5gm
2'-Deoxycytidine hydrochloride Cell Culture Tested $C_9H_{13}N_3O_3.HCl$ MW : 263.68 CAS : 3992-42-5 Assay : $\geq 99\%$ Store below 30°C	TC112-250MG TC112-1G TC112-10G	250mg 1gm 10gm
2-Deoxy-D-ribose Cell Culture Tested $C_5H_{10}O_4$ MW : 134.13 CAS : 533-67-5 Assay : $\geq 99\%$ Store at 2 - 8°C	TC114-1G TC114-5G	1gm 5gm
N6-2'-O-Dibutyryl adenosine-3',5' cyclic monophosphate sodium salt Bucladesine sodium salt, Dibutyryl cAMP sodium salt, Dibutyryl cyclic-AMP sodium salt Cell Culture Tested $C_{18}H_{23}N_5O_8PNa$ MW : 491.37 CAS : 16980-89-5 Assay : $\geq 95\%$ Store at -20°C	TC467-25MG TC467-100MG TC467-250MG TC467-1G	25mg 100mg 250mg 1gm
β-Nicotinamide adenine dinucleotide, From Yeast Cell Culture Tested $C_{21}H_{27}N_7O_{14}P_2$ MW : 663.43 CAS : 53-84-9 Assay : $\geq 95\%$ Store at 2 - 8°C	TC156-100MG TC156-1G TC156-5G TC156-10G	100mg 1gm 5gm 10gm
Thymidine Cell Culture Tested $C_{10}H_{14}N_2O_5$ MW : 242.23 CAS : 50-89-5 Assay : $\geq 99\%$ Store below 30°C	TC174-1G TC174-5G TC174-10G TC174-25G	1gm 5gm 10gm 25gm
Thymine Cell Culture Tested $C_8H_8N_2O_2$ MW : 126.11 CAS : 65-71-4 Assay : $\geq 98\%$ Store below 30°C	TC175-5G TC175-25G	5gm 25gm

Product	Code	Packing
Uracil Cell Culture Tested $C_4H_4N_2O_2$ MW : 112.09 CAS : 66-22-8 Assay : $\geq 99\%$ Store below 30°C	TC178-25G TC178-100G TC178-500G	25gm 100gm 500gm
Uridine Cell Culture Tested $C_9H_{12}N_2O_5$ MW : 244.2 CAS : 58-96-8 Assay : $\geq 99\%$ Store at 2 - 8°C	TC181-5G TC181-25G TC181-50G	5gm 25gm 50gm
Chemicals Meeting Pharmacopeial Testing Specifications		
Acetic acid, Glacial New Meets USP and EP testing specifications CH_3COOH MW : 60.05 CAS : 64-19-7 Store below 30°C	TC574M-1000ML	1000ml
Acetylcysteine N-Acetyl-L-Cysteine New Acetylcysteine (From non-animal source) Meets USP and EP testing specifications $C_6H_9NO_3S$ MW : 163.191 CAS : 616-91-1 Store at 2-8°C	TC054M-5G TC054M-25G TC054M-100G	5gm 25gm 100gm
L-Alanine New (From non-animal source) Meets USP and EP testing specifications $C_3H_7NO_2$ MW : 89.09 CAS : 56-41-7 Store below 30°C	TC056M-25G TC056M-100G TC056M-500G TC056M-1KG	25gm 100gm 500gm 1kg
Ammonium chloride New Meets USP and EP testing specifications ClH_4N MW : 53.49 CAS : 12125-02-9 Store below 30°C	TC092M-100G TC092M-500G	100gm 500gm
L-Arginine New (From non-animal source) Meets USP and EP testing specifications $C_6H_{14}N_4O_2$ MW : 174.2 CAS : 74-79-3 Store below 30°C	TC052M-25G TC052M-100G TC052M-500G TC052M-1KG	25gm 100gm 500gm 1kg
L-Arginine monohydrochloride (From non-animal source) Meets USP and EP testing Specifications $C_6H_{14}N_4O_2.HCl$ MW : 210.66 CAS : 1119-34-2 Store below 30°C	TC053M-25G TC053M-100G TC053M-500G TC053M-1KG	25gm 100gm 500gm 1kg